

Study on the Spectrometer Controlled by Computer Using PIC16F877

Jong Chol Jin, Sok Chol and Ho Tong Gil

The great leader Comrade Kim Il Sung said:

“Therefore, nuclear physicists must intensify their research into atomic energy.”

There is few paper on connecting photo spectrometers with computer for decision of radiation dose intensity using microprocessor PIC16F877.

In this paper, we made a spectrometer for decision of radiation dose connected to computer by using PIC16F877, and considered how to decide absorbance by means of it.

1. The Composition of Circuit of Spectrometer connected to Computer using PIC16F877

In general, the absorbance of any material is proportional to the concentration of that material when transmission length is certain. [1]

$$A = \varepsilon cl$$

where ε is the absorption coefficient in mol, c is density of the material, and l is transmission length.

Based on this equation, the absorbance curve is figured from the values of absorbance of the normal concentration of that material and the concentration of sample is determined from this curve by measuring the sample's absorbance. Diagram of spectrometry “SPECKOL-ZV” is shown as in Fig 1.




Figure1. Diagram of spectrometer

As shown in Figure 1 transmission light is converted to current signal in photoelectric cell. Converting and amplifying this signal into voltage signal, a signal proportional to the absorbance is obtained.

As the earlier devices does not satisfy the linear-

ity of amplitude circuit and the measuring data are displayed on a gauge, measurement error is significant. Furthermore, they should be used in experiments which measures absorbance of according materials, under the condition of known maximum absorption wave of the measured material.

Also the disadvantage is that the experiment process and data processing need lots of effort.

We manufactured the amplifying circuit to measure the signals of wavelength and luminous intensity of the spectrometer at the range of 0~5V and those signals were transformed to digital ones by A-D transforming functionality of PIC.

2. Control device and measurement program.

We implemented a control program, called PIC16F877 for absorbance measured data display and data transmission using USART module according to wave by A-D conversion of measured signals and 7 component LED device by using VC++.

Algorithms of device control and interrupt service are shown in Figure 2.




Figure2. Device control algorithm (a) and interrupt service algorithm (b)

Device control program, firstly, reads value of C port control register PORTC, identifies the current measured values from that and saves A-D converted value of each light intensity read from RA1 pin to corresponding variable if it is set the 0 scales and 100 scales compensation.

Also by substituting the value to Timer0 register, timer interrupt is set and then reading and averaging the ADC value in Timer0 interrupt routine it is calibrated. Main routine of program transmits wavelength and absorbance value continuously with a certain delay.




Figure 3. Measurement algorithm.

In timer0 interrupt service routine, wavelength and light intensity signal are converted to digital data to calculate the average value

Absorbance measuring program of sample using VC++ processes reading data from RS-232 port.

Then interface was coded by using Flash8.0 and linked to VC++, so it is equipped with device connection, data processing, data display functionality with sufficient visual effect and user interface rationally.

Measurement algorithm is shown in Figure 3.

On the interface, the “measure-start” and “measure-end” commands are input from user and send to the main program. To do this, a global value `g_State` is defined and it calls `fs_command()` to send data in the start and end button event handler.

As shown in Figure3, when running the program firstly, on main window it can allow to experiment measuring line generation or measure according to the experiment type.

In the case of unknown samples, data of measurement experiment and its number are input and then the program gets into the measurement state.

Absorbance measurement of sample can be repeated for several times and after terminating the meas-

urement, the corresponding measuring line is selected and the concentration in measurement section item is decided .

With the data browser, the data of measurement line and measurement of sample can be read and here about the experiment data, time, experimenter and content can be browsed in detail.

3. Sample analysis by using photo spectrometer

Using colorimetric analysis of Fe(III) by rhodanide method the reliability of device was confirmed.

This method is based on interaction of Fe_3^+ and SCN^- in acid medium forms red compound.

At this time, the maximum absorption wave of ion Fe(III) was 489nm.

The coloring agent of rhodanide (concentration of SCN^- is 0.623%) was prepared and after drawing the measuring line, the absorbance was measured. (Table1).

In order to confirm characteristics of the device, the halotrichite was prepared in order of about 10^{-12}mol/L .

Table1. Absorbance according to the standard series of Fe(III)

Concentration $/(\cdot 10^{-12}\text{mol}\cdot\text{L}^{-1})$	1	2	3	4	5	6	7
Absorbance	0.105	0.221	0.348	0.437	0.610	0.657	0.725

The change of absorbance Fe(III) decided by obtained measurement value and measuring line is shown in Table2.

Table2. The result of absorbance of Fe samples

No.	1			2		
Absorbance	0.364	0.382	0.321	0.638	0.591	0.684
Concentration Fe(III) $/(\cdot 10^{-12}\text{mol}\cdot\text{L}^{-1})$	3.26	3.42	2.86	5.82	5.38	6.25
Average value $/(\cdot 10^{-12}\text{mol}\cdot\text{L}^{-1})$	3.18 ± 0.3			5.82 ± 0.43		

The relative aberrations of the above results were 9.4% and 7.4%, respectively.

Subsequently, by using above material and measuring by the device connected to computers errors were compared.

The absorbance of Fe(III) ion by spectrometer connected to computer is shown as in Table 3.

Table 3. Analysis results of absorbance of Fe(III) ion by spectrometer connected to computer

No	1			2		
absorbance	0.423	0.418	0.433	0.562	0.571	0.556
concentration /($\cdot 10^{-12} \text{mol} \cdot \text{L}^{-1}$)	3.81	3.76	3.90	5.11	5.19	5.05
Average value /($\cdot 10^{-12} \text{mol} \cdot \text{L}^{-1}$)	3.82 ± 0.07			5.11 ± 0.07		

Relative errors from table 3 are 1.8% and 1.4%, respectively.

Thus, it can be proved that the measurement accuracy can be improved by connecting spectrometer to computer and conduct quantitative analysis of material in range of relative aberration below 2.5%.

Conclusion

1. We manufactured the circuits for the signal of wavelength sensor and the light intensity amplification of the spectrometer. Then we made the device connected to computer.
2. We verified the characteristics of computer- controlled spectrometer through the analysis of Fe(III) ion by rhodanate.
3. The related error in colorimetric analysis of Fe(III) by rhodanide method is less than 2.5%

References

- [1] 조광원 등; 계기 분석, 김일성종합대학출판사, 14~32, 주체96(2007).